



International Journal of Dermatology, Venereology and Leprosy Sciences

E-ISSN: 2664-942X
P-ISSN: 2664-9411
www.dermatologypaper.com
Derma 2021; 4(2): 55-60
Received: 08-07-2021
Accepted: 10-08-2021

Dr. Neha Saraswat
Junior Resident,
Department of D.V.L.
Maharishi Markandeshwar
Medical College & Hospital,
Solan, Himachal Pradesh,
India

Dr. Rajwinder Singh
Professor, Department of
D.V.L. Maharishi
Markandeshwar Medical
College & Hospital, Solan,
Himachal Pradesh, India

Dr. Vinay Shanker
Professor and H.O.D.
Department of D.V.L.
Maharishi Markandeshwar
Medical College & Hospital,
Solan, Himachal Pradesh,
India

Corresponding Author:
Dr. Neha Saraswat
Junior Resident,
Department of D.V.L.
Maharishi Markandeshwar
Medical College & Hospital,
Solan, Himachal Pradesh,
India

Assessment of visual analogue scale (VAS) in patients undergoing elective hair reduction with diode laser

Dr. Neha Saraswat, Dr. Rajwinder Singh and Dr. Vinay Shanker

DOI: <https://doi.org/10.33545/26649411.2021.v4.i2a.94>

Abstract

Background: Diode lasers are the most popular for hair removal due to the deep penetration and targeting of the hair follicle. Our study aimed at improvement in visual analogue scale (VAS), in patients undergoing elective hair reduction with diode laser.

Aim: Assessment of visual analogue scale at every sitting.

Methodology: this prospective observational study was undertaken in the department of dermatology & venereology and leprosy (DVL), Maharishi Markandeshwar Medical College and Hospital, Kumarhatti, Solan during the Dec 2019 to Dec 2020. Total of 30 patients aged more than 16yrs fulfilling inclusion criteria were included in the present study. Lumenis Light Sheer Desire Laser System® of 805nm was used in treatment.

Result: In present study, total of 30 patients fulfilling inclusion criteria were included in the study after obtaining the informed consent. All the patients who participated in study were females. There was a significant higher mean of visual analogue score among the patients compared to the 1st sitting. ($p < 0.05$)

Conclusion: There was a significant improvement in the mean score of VAS among the patients at 5th sitting compared to the 1st sitting. The Diode laser technique is efficient and safe in patients undergoing the elective hair reduction.

Keywords: diode laser and visual analogue scale.

Introduction

Diode lasers are the most popular for hair removal due to the deep penetration and targeting of the hair follicle. Often referred to as the gold standard in hair removal, the diode laser (810nm) is highly reliable, has high papilla absorption, and can address a wide range of skin types, including skin types V and VI.(1) Diode lasers use the principle of selective photolysis to target specific chromophores in the skin, usually melanin or blood. The lasers damage the chromophores by selective heating, without damaging the surrounding tissue. Melanin is the chromophore in the hair follicles, which is targeted and damaged, which results in disrupted growth and regeneration of hair. Our study has evaluated Subject satisfaction with visual analogue scale (VAS) from 0 for complete dissatisfaction to 10 for complete satisfaction.

Material and Method

The present study was undertaken in the department of dermatology & venerology and leprosy (DVL), Maharishi Markandeshwar Medical College and Hospital, Kumarhatti, Solan.

Study Design: Prospective Observational Study.

Duration: December 2019 to December 2020

Selection of patients

All the patients, who came to the Out Patient Department of Dermatology, Venereology and leprosy of MMMC&H, for the first time, for elective hair reduction with diode laser.

Sample size: Minimum sample size of 30.

Inclusion criteria

Age > 16yrs

Exclusion criteria

- History of previous laser hair removal treatment
- Current or past history of skin cancer
- Severe concurrent medical conditions
- Any active cutaneous eruption in the treatment area

Procedures:

Preparation:

1. Detailed clinical history to rule out infections, drug intake, etc.
2. In patients of hirsutism modified Farriman-Gallwey score was calculated, along with following investigations to rule out the cause and treated accordingly
 - Serum testosterone
 - Serum prolactin
 - USG pelvis
 - Serum T3, T4 and TSH
3. Avoid hair removal by plucking, bleaching or waxing for 6 weeks before the procedure.
4. Sunscreens were prescribed 4 to 6 weeks before the procedure, especially on exposed areas.
5. Count number of hairs in per cm² area
6. Shave the area before the procedure.
7. Obtain an informed consent.
8. Peri-follicular oedema and erythema is an expected side effect, and this was explained to the patient.

Equipment

Lumenis Light Sheer Desire Laser System® OF 805nm, used in Dermatology, Venereology and Leprosy department of MMMCH, Kumarhatti, Solan.

Technique:

1. Clean the area to remove oils from the skin
2. Position the patient comfortably.
3. Use eye protection for both the operator as well as the patient
4. Apply ice pack on area to be treated
5. Place the handpiece (XC or HS tip), perpendicular to the skin according to area to be treated. Ensure that it is pressed down gently and firmly to bring the follicles close to the surface and temporarily decrease the blood supply to the local area.
6. A total of 5 sittings were given at an interval of 4 to 6 weeks.

Post procedure care

1. Although most machines have inbuilt cooling systems, ice packs or cooling agents was used post-procedure to minimize pain and swelling.
2. In case of inadvertent epidermal injury, topical antibiotics was used.

3. Strict sun protection and the use of sunscreens.
4. Topical steroid creams were prescribed, if there were any side effects like erythema and swelling.

Pulse width and fluence were used as per the manufacturer’s guidelines for skin type IV and V For XC tip:

Sitting	Pulse Width (Ms)	Fluence (J/Cm2)
First	100	18-20
Second	100	20-22
Third	100	22-24
Fourth	30	20-22
Fifth	30	22-24

For HS tip

Sitting	Pulse width (ms)	Fluence (j/cm2)
first	100	5.0-5.5
second	100	5.5-6.0
third	30	5.0-5.5
fourth	30	5.5-6.0
fifth	30	6.0-6.5

Ethical consideration: Informed and written consent (in the language they best understand) was taken from each subject before collecting data and blood sample. Only those individuals, who volunteer to participate in the study, were included and the data was kept confidential. The study did not impose any burden on the subjects and the Institute; therefore, the study is ethically justified.

Results

In present study, total of 30 patients fulfilling inclusion criteria were included in the study after obtaining the informed consent. All the patients who participated were females, as males are not as concerned as females about their hairs. The mean age of the patients was found to be 27.83±6.53yrs of age, with majority of patients in group of 17-25yrs (40%) followed by 26-30yrs of age (33.3%).

Table 1: Mean age of the patients in present study

	N	Minimum	Maximum	Mean	SD
Age in yrs	30	17	43	27.66	6.47

Out of total 30 patients, minimum age of the patient was 17yrs and maximum age was 43yrs. Mean age of the patients was found to be 27.66±6.47yrs

Table 2: Showing the age wise distribution of patients

	Frequency	Percent
17-25yrs	12	40
26-30yrs	10	33.3
>31yrs	8	26.7

Majority of patients in group of 17-25yrs (40%), followed by 26-30yrs of age (33.3%) and >30yrs in 26.7%.




Fig 1: Showing the age wise distribution of patients
 Out of 30 patients, 12 patients were between the age of 17-25yrs, 10 patients between 26-30yrs and 8 patients were of age >31yrs.

Table 3: Showing the mean level of visual analogue scale score during various sitting

	N	Minimum	Maximum	Mean	SD
1st sitting Visual analogue scale	30	.0	.0	.000	.00
2nd sitting Visual analogue scale	30	1.0	5.0	2.00	1.01
3rd sitting Visual analogue scale	30	1.0	6.0	3.70	1.34
4th sitting Visual analogue scale	30	2.0	8.0	5.46	1.35
5th sitting Visual analogue scale	30	2.0	9.0	6.36	1.51

Subject satisfaction was assessed using visual analogue scale (VAS) from 0 for complete dissatisfaction to 10 for complete satisfaction. Mean VAS at 1st sitting; 0, 2nd sitting; 2.0, 3rd sitting; 3.70, 4th sitting; 5.46, 5th sitting; 6.36.




Fig 2: Showing the mean level of visual analogue scale score
 VAS increasing from 0 at 1st sitting to 6.36 at 5th sitting

Table 4: Showing the mean level of hair density during various sitting

	N	Minimum	Maximum	Mean	SD
1st sitting Hair Density	30	3.0	44.0	17.33	10.41
2nd sitting Hair Density	30	3.0	37.0	13.86	8.75
3rd sitting Hair Density	30	1.0	30.0	11.06	7.74
4th sitting Hair Density	30	1.0	22.0	8.26	5.75
5th sitting Hair Density	30	1.0	20.0	6.90	4.97

Hair density was measured by drawing 1x1 cm² box and calculating the terminal hairs in that box. Mean hair density at 1st sitting; 17.33, 2nd sitting; 13.86, 3rd sitting; 11.06, 4th sitting; 8.26, 5th sitting; 6.90.




Fig 3: Showing the mean level of hair density
Mean hair density decreasing from 17.33 at 1st sitting to 6.90 at 5th sitting

Table 5: Comparison of mean visual analogue scale score at various sittings using paired t-test

Paired Samples Statistics				
		Mean	SD	Paired t-test (p-value)
Pair 1	1st sitting Visual analogue scale	.00	.00	<0.0001
	2nd sitting Visual analogue scale	2.00	1.02	
Pair 2	1st sitting Visual analogue scale	.00	.00	<0.0001
	3rd sitting Visual analogue scale	3.70	1.34	
Pair 3	1st sitting Visual analogue scale	.00	.00	<0.0001
	4th sitting Visual analogue scale	5.47	1.36	
Pair 4	1st sitting Visual analogue scale	.00	.00	<0.0001
	5th sitting Visual analogue scale	6.37	1.52	

P value was highly significant with every subsequent sitting on comparing with first sitting.




Fig 4: Comparison of mean visual analogue scale score at various sittings

Table 6: Comparison of mean hair density at various sittings using paired t-test

		Mean	SD	Paired t-test (p-value)
Pair 1	1st sitting Hair Density	17.33	10.42	<0.0001
	2nd sitting Hair Density	13.87	8.75	
Pair 2	1st sitting Hair Density	17.33	10.42	<0.0001
	3rd sitting Hair Density	11.07	7.75	
Pair 3	1st sitting Hair Density	17.33	10.42	<0.0001
	4th sitting Hair Density	8.27	5.75	
Pair 4	1st sitting Hair Density	17.33	10.42	<0.0001
	5th sitting Hair Density	6.90	4.98	

P value was highly significant with every subsequent sitting on comparing with first sitting.




Fig 5: Comparison of mean hair density at various sittings

Discussion

Unwanted facial hair (UFH) is an important issue: over 40% of the women in the general population have some degree of UFH, ranging from nearly invisible growth to extensive growth. It has psychological and psychosocial impact, independent from the severity, which should not be underestimated. Hirsutism is the most common cause of unwanted facial hairs, which is defined as the presence of excessive terminal hair in females in androgen dependent areas.

The diode laser system generates a laser light of ~800 nm wavelength. In a prospective controlled trial, where they used this laser, resulted in significant reduction of hair growth. Overall, clinical studies with the diode laser for hair reduction have reported variable success rates, ranging from 65-75% hair reduction at 3 months after one to two treatments with fluences of 10-40 J/cm² to >75% hair reduction in 91% of subjects 8 months after three to four treatments at 40J/cm².

Alternatively, conflicting data exists that, whether the fluence of the laser, should be between 15-40J/cm² for therapeutic success. In darker skin, less fluence may be required to attain the same efficacy. In general, the diode lasers are found to be better tolerated by people of darker skin types (V-VI) than the ruby laser; because the energy absorbed by the melanin is ~30% less than that of ruby laser, but has a better optical penetration. Although the incidence of pigmentary skin changes has been as high as 29%, most of these changes are transient and well tolerated.

Advantages of diode laser are as follows:

- The Diode laser features a longer wavelength and thus able to provide better results among dark-skinned people.
- This laser provides safer and deeper penetration into the skin layer.
- Large areas of the body tend to recover faster.

Treatment with diode laser is effective on medium to dark hair on any skin type, when the hair is in its anagen phase of

growth; because of this reason, multiple sessions at an interval of one month are required for the treatment, as at a time not all the hairs are in anagen phase of growth. Therefore; the number of sessions required for the hair reduction varies.

The mean age of patients in the present study was 27.66±6.47yrs. Similar to present study, Baig T *et al.*, also documented the mean age of participants to be 26.12±5.83 years. (2). In an another study women were between the ages 24 and 41 years with a mean age of 32.27±4.21 years, which is slightly higher than our study.(3)

VAS is a subject satisfaction score; with each sitting, there was a significant higher mean of visual analogue score among the patients compared to the 1st sitting. (*p*<0.05) The mean VAS score increased from 2.00±1.01 at 2nd sitting, 3.70±1.34 at 3rd sitting, 5.46±1.35 at 4th sitting and 6.36±1.51 at 5th sitting. Similar to our study ayatollahi *et al* observed increase in VAS from 0 to 5.07 after diode laser hair reduction [4].

With the treatment, hair density was found to reduce significantly at each sitting compared to the mean of 1st sitting. (*p*<0.05) The mean level reduced from 17.33±10.41 at 1st sitting, 13.86±8.75 at 2nd sitting, 11.06±7.74 at 3rd sitting, 8.26±5.75 at 4th sitting and 6.90±4.97 at 5th sitting. In a study conducted over 14 female patients, calculated hair density within 2x2 cm² area. Average hair density at baseline was 69.1, which came down to 20.6 at 6 months [5]. HR-AL hamamy observed hair density to be 21.58±11.59 at baseline, 15.04±9.17 after 1st session, 11.77±7.55 after 2nd session, 8.96±6.77 after 3rd session [6]. Similarly, in a study over 8 patients, approximate decrease in hair growth was from 75%-90% [7]. In a study conducted over people of various ethnicity, per cm² reduction of hair was noticed, 92.5% for Asians, 81.5% for people with black skin,86.9% for mixed race and 88.2% for people of white coloured skin [8]. Zhou *et al* calculated hair density by drawing a circle of 3.5mm and calculating the hairs within the circle. Reduction rate of hair growth noticed was 72±18 at the end of 3

months^[3].

Most of the patients had erythema and peri-follicular edema, which was transient. 1st degree burn was noticed in one (0.03%) patient at chin on 2nd sitting, and was given topical anti-inflammatory.

Conclusion

This prospective observational clinical study of 30 female patients confirms that diode laser combined with skin cooling provides a safe and effective method for long-term hair reduction in Fitzpatrick skin types IV and V. By using conservative parameters, correct hormonal workup and ultrasonography, careful patient selection and adequate gap between two sessions, gave us an excellent result with no long-term side effects and good patient satisfaction rates. To conclude, the Diode laser technique is efficient and safe in patients undergoing the elective hair reduction. But various multi-centric studies with a follow-up of more than 1 year are required to determine the long term effects of diode laser hair removal.

References

1. Courtney E, Goldberg D. Clinical Evaluation of Hair Removal Using an 810 nm Diode Laser with a Novel Scanning Device. *J Drugs Dermatol*. 2016;15(11):1330-3.
2. Ekbäck MP, Lindberg M, Benzein E, Årestedt K. Social support: An important factor for quality of life in women with hirsutism. *Health Qual Life Outcomes* 2014;12(1):183.
3. Gaber M, El-Sayed S. Quality of life in patients having hirsutism. *Menoufia Med J* 2021;34(2):477-81.
4. Ayatollahi A, Samadi A, Rajabi-Estarabadi A, Yadangi S, Nouri K, Firooz A. Comparison of efficacy and safety of a novel 755-nm diode laser with conventional 755-nm alexandrite laser in reduction of axillary hairs. *Lasers in medical science*. 2020;35(2):373-8.
5. JD DJ. Clinical evaluation of hair removal using an 810 nm diode laser with a novel scanning device. *Journal of drugs in dermatology*. 2016;15(11):1330-3.
6. AL-Hamamy HR, Saleh AZ, Rashed ZA. Evaluation of Effectiveness of Diode Laser System (808 nm) versus Intense Pulse Light (IPL) in the Management of Unwanted Hair: A Split Face Comparative Study. *Int J Med Physics, Clin Eng Radiat Oncol* 2015;04(01):41-8.
7. Greppi I. Diode laser hair removal of the black patient. *Lasers in Surgery and Medicine: The Official Journal of the American Society for Laser Medicine and Surgery* 2001;28(2):150-5.
8. Zaleska I, Atta-Motte M. The Effectiveness of Diode Laser 805 nm Hair Removal in Groups of Various Ethnicity. *J Med-Clin Res & Rev* 2018;2(4):1-9